


Upgrade path to four position for your Smart Evaporator C1

4Kit for Smart Evaporator C1

“4kit” transforms your C1 into FOUR position simply by exchanging the head parts and heat block.



Smart Evaporator C 1

Smart Evaporator C1 with 4kit attached


4kit is a set of:

- 1 head part that accommodates four plugs (to be connected directly to the C1 plug adaptor)
- 4 each of P1 and P2 Spiral Plugs
- 4 sealing plugs
- 1 heat block that replaces the aluminum beads (one block can accommodate any of below vials)
- 4 each of 2ml, 4ml and 6ml vials



Four plugs



Heat block



Please check the type of your C1 in order to choose the fitting 4kit.

Smart Evaporator C1	4Kit	
 <p>Type of attachable Smart Evaporator: Old Model (2012.3~2018.3)</p>	Part number	4kit UW4-C1
	Heat block size	φ 50mm
	Compatible containers	2ml ,4ml, 6ml vials
	Compatible samples	4channels
 <p>Type of attachable Smart Evaporator: New Model (2018.4~)</p>	Part number	4kit UW4-C1B
	Heat block size	φ 69mm
	Compatible containers	2ml ,4ml, 6ml vials
	Compatible samples	4channels

■ Distributor

■ Manufacturer

BC BioChromato

BioChromato, Inc.

1-12-19 Honcho, Fujisawa, Kanagawa-ken 251-0053JAPAN
 TEL: +81-466-23-8382 Fax: +81-466-23-8279
 https://www.bicr.biz/ E-mail: info@bicr.co.jp