

RAPID Tube Plate 600

INSTRUCTION FOR HOW TO USE RAPID Tube Plate 600

Thank you for purchasing RAPID Tube Plate 600.
Please read this manual thoroughly before using the
RAPID Tube Plate 600. Please keep this manual
somewhere easily accessible.



Product Specification

Catalog Code	RTP6001-20/RTP6001-160	
Product Name	RAPID Tube Plate 600	
Material	Polypropylene	
Volume	600 μ L/well	
Functional Temperature Range	-80°C~130°C	
Heat Deflection Temperature	90°C~130°C	
Density	0.9 mg/m ²	
Water Absorbency	<0.01%	
Units	RTP6001-20 (20 plates)/RTP6001-160 (160plates)	
No Sterile	Sterilizable in Autoclave (121°C, 20 min, 1 atm)	Human DNA, DNase/RNase free

※ Contents are subject to change without prior notice.

Caution

- For research use only
- Observe the maximum filling volume.
- Take appropriate safety measures to prevent the escape of substances in the plates during storage and transportation.
- Only use undamaged and visibly perfect well plate from intact packages.
- Store at room temperature and protect from direct sunlight and UV light.
- Once a package is opened, we do not guarantee the quality of the product.

Manufacturer



BioChromato, Inc.

1-12-9 Honcho, Fujisawa, Kanagawa-ken 251-0053 JAPAN

TEL : +81-466-23-8382 FAX : +81-466-23-8279

E-mail : info@bier.co.jp

Website: <https://www.biochromato.com>