


# Do you have sample vaporization issues when culturing cells in 96-well plate?

## RAPID Slit Seal

A seal that can prevent volatilization even when punched with an autosampler!



### Features

- The seal (slit) returns to its original state after insertion and prevents solvent volatilization
- No adhesive above wells so sample leaching is avoided (contamination prevented)
- **Sterilized version** (applicable to cell culture purposes)

## Sales of STERILIZED version started

We are pleased to announce that RAPID Slit Seal is now available in GAMMA - STERILIZED version.

**Useful for edge effect reduction, e.g. when culturing in 96-well plate!**

Product Name:	Sterilized RAPID Slit Seal (10 sheets)
Quantity:	10 sheets / box
Model Number:	RSS-S96-10-S
Price:	Please ask your distributor.


Product Name:	Sterilized RAPID Slit Seal (50 sheets)
Quantity:	50 sheets (10sheets ×5 boxes)
Model Number:	RSS-S96-50-S
Price:	Please ask your distributor.

**RAPID Slit Seal**

検索



FAQ · Evaluation Data



Evaporation rate result is almost the same with the one of the general plate seal.

Size : W 80mm×L 122mm

Material : PET, Silicon

Operating temperature : -80°C~100°C\* \*Based on the evaluation test carried out at 100°C for 30 minutes.

■ Distributor

■ Manufacturer

**BC BioChromato**

**BioChromato, Inc.**

1-12-19Honcho, Fujisawa, Kanagawa-ken 251-0053 JAPAN

TEL: +81-466-23-8382 Fax: +81-466-23-8279

https://www.bicr.biz/

E-mail: info@bicr.co.jp